

Disinfection and cleaning plan including skin protection plan for dental surgeries

in accordance with DGUV Regulation 1 (German Statutory Accident Insurance Association), TRBA 250 and Section 36 of the German Infectious Diseases Protection Act



What? The object to be maintained	How? Type of maintenance	With what? Work materials (e.g. disinfectants, VAH ¹ -listed), preparation, dosage, reaction time ¹⁾	When? Timing, frequency, sequence of hygiene measures	Who? Name and signature of the person responsible																					
Instruments																									
<p>Instruments and aids for</p> <ul style="list-style-type: none"> general, preventive, restorative or orthodontic treatments, such as impression trays, trays, cassettes, tablets, glasses, suction cannulas or detachable suction couplings surgical, periodontal or endodontic treatments, e.g. syringes, forceps, levers <p>Rotating or oscillating instruments for</p> <ul style="list-style-type: none"> general, preventive, restorative or orthodontic treatments, e.g. polishers, drills, milling tools, grinders, finishing burs surgical, periodontal or endodontic treatments, e.g. endodontic/periodontal instruments and surgically rotating instruments 	<p>Instrument reprocessing:</p> <ol style="list-style-type: none"> Proper preparation of the reprocessing and contamination-proof transport to the reprocessing location Removal of coarse organic soiling with cellulose Cleaning and disinfection <p>3.1. Wet-chemical:</p> <ul style="list-style-type: none"> Immersion of the instruments (disassembled or with joints opened), bubble-free, covered inside and out, in a disinfecting and cleaning solution, mechanical cleaning (ultrasonic bath if needed) intermediate rinsing with water Chemical disinfection, rinsing after the end of the reaction time with suitable water, drying <p>3.2. Machine:</p> <ul style="list-style-type: none"> In a washer-disinfector: disinfect, clean and dry in accordance with the manufacturer's information <ol style="list-style-type: none"> Check for cleanliness, integrity, if needed repeat step 3 Identification Care and maintenance Function test Labelling (as required) Steam sterilisation with an effective and validated procedure: <ol style="list-style-type: none"> Non-critical instruments^{**}; steam sterilisation not required Semi-critical instruments^{**}; unpackaged on trays, sieve trays, tablets, stands, bowls, etc. or packed in sterile product packaging in accordance with EN 868 Critical instruments^{***}; packaged in sterile product packaging in accordance with EN 868 Documented release of the medical devices for use Contamination-proof storage, storage time in accordance with EN 868 	<table border="1"> <thead> <tr> <th></th> <th>acc. to VAH¹⁾</th> <th>tuberculocidal</th> <th>in the ultrasonic bath bacteri-cidal/yeasticidal</th> <th>in the ultrasonic bath incl. tuberculocidal</th> </tr> </thead> <tbody> <tr> <td>ID 212</td> <td rowspan="3">2% / 5 min</td> <td rowspan="3">2% / 60 min</td> <td>2% / 5 min.</td> <td>2% / 60 min</td> </tr> <tr> <td>ID 212 forte</td> <td>2% / 5 min.</td> <td>4% / 30 min</td> </tr> <tr> <td>ID 213</td> <td>2% / 60 min</td> <td>2% / 30 min</td> </tr> <tr> <td>ID 220</td> <td>100% / 1 min</td> <td></td> <td colspan="2">100% / 0,5 min</td> </tr> </tbody> </table> <p>For general and surgical instruments: Disinfection: ID 212, ID 212 forte, ID 213, 2% (20 ml/l water) Cleaning: ID 215²⁾, 2% (20 ml/l water), 15 min.</p> <p>For rotating instruments: Disinfection: ID 220, undiluted Cleaning: ID 215²⁾, 2% (20 ml/l water), 15 min.</p> <p>Washer-disinfector: Program: Temperature/hold time:</p> <p>Steam steriliser: Hygoclave 50 Program: Temperature/hold time:</p>		acc. to VAH ¹⁾	tuberculocidal	in the ultrasonic bath bacteri-cidal/yeasticidal	in the ultrasonic bath incl. tuberculocidal	ID 212	2% / 5 min	2% / 60 min	2% / 5 min.	2% / 60 min	ID 212 forte	2% / 5 min.	4% / 30 min	ID 213	2% / 60 min	2% / 30 min	ID 220	100% / 1 min		100% / 0,5 min		<p>Processing is performed directly after the treatment, taking into account the risk assessment and classification of the medical devices</p> <p>dry storage in washer-disinfector before the disinfection, not over 6 hrs.</p>	
	acc. to VAH ¹⁾	tuberculocidal	in the ultrasonic bath bacteri-cidal/yeasticidal	in the ultrasonic bath incl. tuberculocidal																					
ID 212	2% / 5 min	2% / 60 min	2% / 5 min.	2% / 60 min																					
ID 212 forte			2% / 5 min.	4% / 30 min																					
ID 213			2% / 60 min	2% / 30 min																					
ID 220	100% / 1 min		100% / 0,5 min																						
<p>Transmission instruments for general, preventive, restorative, orthodontic, surgical, periodontal and endodontic treatments, e.g. turbine, hand and angle pieces</p> <p>Auxiliary devices with/without leakage of fluids and/or air or particles, e.g. dental clearance equipment, powder jet equipment, intraoral camera, etc.</p>	<p>Processing is performed in the same way as for the instruments described above (steps 3 and 9 according to the manufacturer information)</p> <p>External: wipe disinfection and cleaning</p> <p>Interior: cleaning Care (lubrication of interior parts)</p> <p>Automated processing: Disinfection/cleaning e.g. in a washer-disinfector Internal cleaning and maintenance of internal parts</p>	<p>FD 322 premium wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333 wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333 forte premium wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 350, FD 350 green disinfecting wipes, undiluted, 1 min. in accordance with VAH⁴⁾</p> <p>Unit: Program: Temperature/hold time:</p>	<p>After every treatment and after contamination</p>																						
Surfaces																									
<p>Surfaces and floors in examination and treatment rooms (e.g. treatment unit, devices, cabinets)</p>	<p>Disinfect and clean surfaces with a mopping and wiping disinfection procedure (using an appropriate wet wiping method). No rinsing or wiping dry before the end of the reaction time. Re-use the surfaces when they are visibly dry.</p>	<p>FD 300, 1% (10 ml/l water), 15 min. acc. to VAH⁴⁾ FD 312²⁾, 1% (10 ml/l water), 15 min. acc. to VAH⁴⁾ FD 312 wet wipes³⁾, undiluted, 15 min. acc. to VAH⁴⁾ FD multi wipes, FD multi wipes compact, FD multi wipes green, FD multi wipes compact green</p>	<p>Every working day (e.g. at the end of the treatment day), additionally as required</p>	<p>Cleaner</p>																					
<p>Objects and surfaces near the patient that have been contaminated and cannot be disinfected thermally or by immersion in solutions (e.g. treatment unit, devices, work surfaces, equipment parts, X-ray equipment, lamp handles, intraoral film shells)</p>	<p>Wiping with quick disinfectants Limit spray disinfection to areas that cannot be reached by wipe disinfection.</p> <p>Ensure complete wetting, comply with reaction time, allow to dry</p> <p>Cleaning and care (e.g. treatment chair)</p>	<p>FD multi wipes, FD multi wipes compact, FD multi wipes green, FD multi wipes compact green FD 322, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333 forte, undiluted, 1 min. in accordance with VAH⁴⁾ FD 366 sensitive, undiluted, 1 min. in accordance with VAH⁴⁾ FD 322 premium wipes, FD 322 top wipes, FD 322 wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333 wipes, FD 333 top wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333 forte premium wipes, FD 333 forte wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 350, FD 350 green disinfecting wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 366 sensitive top wipes, FD 366 sensitive wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 360 cleaning for synthetic leathers, undiluted.</p>	<p>In the treatment area after every treatment and after contamination</p> <p>As required</p>																						
Skin protection plan																									
<p>Requirements are an easily accessible hand-washing station with warm and cold running water, direct dispensers with washing lotion that is gentle on the skin, hand disinfectant and appropriate skin care products as well as single-use towels. No jewellery, watches or rings may be worn on the hands or forearms.</p>																									
<p>Skin, when affected by</p> <ul style="list-style-type: none"> water-soluble working materials, such as disinfectants, cleaning agents, X-ray development chemicals non water-soluble materials, such as organic solvents, petroleum ether accumulation of moisture and softening of skin (maceration), e.g. when wearing gloves 	<p>Cleaning hands: Under running water with gentle liquid soap from a dispenser, dry with single-use towels</p> <p>Skin protection and care: Apply carefully between fingers, paying special attention to fingernail folds and the back of the hands</p>	<p>Gentle liquid washing preparation: HD 435 single-use towel from direct dispenser</p> <p>Protection emulsion W/O Preparation:</p>	<p>Before starting work, after finishing work, in case of visible soiling, before meals, following every visit to the toilet, after sneezing, blowing one's nose or coughing</p>	<p>All employees</p>																					
<p>Skin, in the event of non-physiological colonisation with e.g. bacteria, fungi, viruses</p>	<p>Hygienic hand disinfection: Hand disinfectant in accordance with the rub-in process CEN/EN 1500 to be rubbed in until dry. Hands must be air-dry before putting on the medical single-use gloves.</p> <p>Surgical hand disinfection: Wash hands and forearms to the elbow for 1 minute, brushing fingernails and nail folds if required. Moisten hands and forearms completely by rubbing in the hand disinfectant solution in portions (with the rub-in method as per EN 12791) and keep them moist during the reaction time. Hands must be air-dry before putting on sterile medical single-use gloves.</p>	<p>HD 410, at least 3 ml, 15 sec. acc. to VAH⁴⁾</p> <p>HD 410, at least 10 ml, 1.5 min. acc. to VAH⁴⁾</p>	<p>Before preparing for work, before and after each treatment, after interruptions in treatment, after workstation maintenance, before putting on or taking off gloves</p> <p>Before extensive dental-surgical or oral-surgical procedures, before all dental-surgical or oral-surgical procedures on patients with an increased risk of infection</p>	<p>All employees</p>																					
Special areas																									
<p>Suction unit, suction hoses (internal)</p> <p>Spittoon (internal)</p> <p>Filter</p> <p>Amalgam separator</p>	<p>Cleaning by aspirating a water-air mix, especially when blood has been aspirated</p> <p>1) Disinfection: slow suction of a mix of air and disinfection/cleaning solution</p> <p>2) Cleaning: removal of deposits (calcium, bacterial contamination, salts from jet-stream equipment)</p> <p>3) Disinfection and cleaning of the spittoon</p> <p>Filter change or cleaning in accordance with the manufacturer's specifications</p> <p>Replacement/drainage of the collecting container and disposal of the separated material in accordance with the manufacturer's information</p>	<p>Cold water (at least 1/2 litre)</p> <p>Cold water (approx. 2 litres) concentrate in the care system OroCup Mix and suction off 1) Orotal[®] plus, 2% (20 ml/l water), 5 min. acc. to VAH⁴⁾ Orotal[®] plus pH 7, 2% (20 ml/l water), 30 min. acc. to VAH⁴⁾ Orotal[®] ultra, 1% (10 g/l water), 5 min. acc. to VAH⁴⁾</p> <table border="1"> <thead> <tr> <th></th> <th>Concentration</th> <th>Contact time</th> </tr> </thead> <tbody> <tr> <td rowspan="2">2a)</td> <td>MD 555 cleaner</td> <td>5% (50 ml/l water) 30-120 min</td> </tr> <tr> <td>MD 555 cleaner organic</td> <td>5% (50 ml/l water) 30-120 min</td> </tr> <tr> <td rowspan="2">2b)</td> <td>MD 555 cleaner</td> <td>5% (50 ml/l water) 5-15 min</td> </tr> <tr> <td>MD 555 cleaner organic</td> <td>5% (50 ml/l water) 5-15 min</td> </tr> </tbody> </table> <p>3) MD 550, undiluted</p> <p>Use PPE (personal protective equipment)</p> <p>Do not touch contaminated pieces of the separator; use PPE</p>		Concentration	Contact time	2a)	MD 555 cleaner	5% (50 ml/l water) 30-120 min	MD 555 cleaner organic	5% (50 ml/l water) 30-120 min	2b)	MD 555 cleaner	5% (50 ml/l water) 5-15 min	MD 555 cleaner organic	5% (50 ml/l water) 5-15 min	<p>After every treatment</p> <p>before each disinfection or cleaning</p> <p>1) 1 to 2 times per day, minimum at the end of the working day</p> <p>2a) at least 2 times per week, prior to the lunch break</p> <p>2b) 5 times weekly, preferably before the start of treatment</p> <p>3) Additionally as required</p> <p>Use PPE As required</p> <p>As required</p>									
	Concentration	Contact time																							
2a)	MD 555 cleaner	5% (50 ml/l water) 30-120 min																							
	MD 555 cleaner organic	5% (50 ml/l water) 30-120 min																							
2b)	MD 555 cleaner	5% (50 ml/l water) 5-15 min																							
	MD 555 cleaner organic	5% (50 ml/l water) 5-15 min																							
<p>Suction hoses (external) etc.</p> <p>Spittoon (internal and external)</p>	<p>Disinfection and cleaning through wiping</p>	<p>FD multi wipes, FD multi wipes compact, FD multi wipes green, FD multi wipes compact green FD 322, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333 forte, undiluted, 1 min. in accordance with VAH⁴⁾ FD 366 sensitive, undiluted, 1 min. in accordance with VAH⁴⁾ FD 322 premium wipes, FD 322 top wipes, FD 322 wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333 wipes, FD 333 top wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 333 forte premium wipes, FD 333 forte wipes, undiluted, 1 min. in accordance with VAH⁴⁾ FD 350, FD 350 green disinfecting wipes, undiluted, 1 min. acc. to VAH⁴⁾ FD 366 sensitive top wipes, FD 366 sensitive wipes, undiluted, 1 min. in accordance with VAH⁴⁾</p>	<p>After every treatment</p>																						
<p>Impressions e.g. alginate, silicones, polyether rubber, hydrocolloids, polysulphides, etc. as well as workpieces and equipment e.g. prostheses, bridges, etc., plaster models, bite registrations and bite plates made of wax, shellac or plastics, as well as articulators, occluders, fixators, face bows and more.</p>	<p>Disinfection and cleaning with Automated procedures: in accordance with the manufacturer's specifications</p> <p>Wet-chemical procedures: Disinfection and cleaning by soaking in solution, after the advised reaction time rise under running water and air dry if necessary</p>	<p>Device/manufacturer: Hygojet Preparation: MD 520, undiluted, 5 min. acc. to VAH⁴⁾</p> <p>MD 520, undiluted, 5 min. acc. to VAH⁴⁾</p>	<p>After removal from mouth prior to transport to dental laboratory after receiving from dental laboratory</p>																						
<p>Objects with acid-soluble coatings e.g. cement, calculus, oxide film, flux residues</p>	<p>Cleaning (e.g. in an ultrasonic bath)</p>	<p>MD 530, undiluted</p>	<p>As required</p>																						
<p>Objects with coatings of plaster or alginate</p>	<p>Cleaning (e.g. in an ultrasonic bath)</p>	<p>MD 535, undiluted</p>	<p>As required</p>																						
Miscellaneous																									
<p>Oral antiseptics</p>	<p>Rinse the oral cavity or soak swabs with antiseptic to reduce germs in the oral cavity of the patient and to safeguard the surgery team against infection at the same time</p>	<p>Preparation, dosage: Lunos[®] Dental rinse, undiluted (20 ml), 1 min.</p>	<p>Before surgical procedures, for patients with an increased risk of infection</p>	<p>Patient</p>																					
<p>Practice laundry Textile practice clothing, professional clothing, towels, wiping cloths, etc.</p>	<p>Practice laundry is to be collected in sufficiently resistant and thick containers/sacks and reprocessed with thermal or chemical-thermal wash and disinfection procedures</p>	<p>Device/manufacturer: Wash program/temperature: Preparation/dosage:</p>	<p>After use After contamination</p>	<p>Operator</p>																					
<p>Waste Waste from treatment and examination rooms</p>	<p>Pointed, sharp, fragile objects, body parts (e.g. teeth, tissue) with blood, secretions from contaminated refuse (e.g. saliva ejector, swabs, etc.), as well as waste contaminated with agents of notifiable diseases must be immediately collected at their point of origin in non-penetrable, tear- and moisture-proof and sealed containers and must be transferred to secure sealed containers for burning. Lead foil, old films, fixing baths, developer, residues containing mercury, old medications and batteries are disposed of by specialist disposal firms. Household-like waste is collected in waste bags and sorted and disposed of in accordance with the local guidelines.</p>																								

^{*} Non-critical medical devices: medical devices that merely come into contact with intact skin
^{**} Semi-critical medical devices: medical devices that come into contact with mucosa or pathologically affected skin
^{***} Critical medical devices: medical devices that penetrate skin or mucosa and come into contact with blood, internal tissue or organs (including wounds) in the process, and medical devices for application of blood, blood products and sterile drugs. Patients suffering from or thought to suffer from clinical conditions of transmissible spongiform encephalopathy (CJK, vCJK etc.) are to be treated in institutions that are sufficiently equipped to prevent infection. The reprocessing of reusable instruments and materials must only be performed in institutions with an externally certified quality management system in accordance with DIN EN ISO 13485/13488.

¹⁾ For further information, refer to the product information leaflets available on the Internet.
²⁾ In accordance with EN 13727, EN 13624, EN 14561, EN 14562
³⁾ Use disinfectants safely. Always read the label and product information before use.
⁴⁾ VAH-listed: bactericidal/yeasticidal

This disinfection and cleaning plan was provided by
 DÜRR DENTAL SE · Höpfigheimer Str. 17
 74321 Bietigheim-Bissingen · Germany · www.duerrdental.com
 Tel. +49 71 42/7050 · Fax +49 71 42/705500



Stamp of the surgery/date/signature of the operator

Disinfection and cleaning plan including skin protection plan for dental surgeries

in accordance with DGUV Regulation 1 (German Statutory Accident Insurance Association), TRBA 250 and Section 36 of the German Infectious Diseases Protection Act



What? The object to be maintained	How? Type of maintenance		When? Timing, frequency, sequence of hygiene measures	Who? Name and signature of the person responsible
Instruments				
<p>Instruments and aids for</p> <ul style="list-style-type: none"> general, preventive, restorative or orthodontic treatments, such as impression trays, trays, cassettes, tablets, glasses, suction cannulas or detachable suction couplings surgical, periodontal or endodontic treatments, e.g. syringes, forceps, levers <p>Rotating or oscillating instruments for</p> <ul style="list-style-type: none"> general, preventive, restorative or orthodontic treatments, e.g. polishers, drills, milling tools, grinders, finishing burs surgical, periodontal or endodontic treatments, e.g. endodontic/periodontal instruments and surgically rotating instruments 	<p>Instrument reprocessing:</p> <ol style="list-style-type: none"> Proper preparation of the reprocessing and contamination-proof transport to the processing location Removal of coarse organic soiling with cellulose Cleaning and disinfection Wet-chemical: <ul style="list-style-type: none"> Immersion of the instruments (disassembled or with joints opened), bubble-free, covered inside and out, in a disinfecting and cleaning solution, mechanical cleaning (ultrasonic bath, if needed) intermediate rinsing with water Chemical disinfection, rinsing after the end of the reaction time with suitable water, drying Mechanical: <ul style="list-style-type: none"> In a washer-disinfector: disinfect, clean and dry in accordance with the manufacturer's information Check for cleanliness, integrity, if needed repeat step 3 Identification Care and maintenance Function test Labelling (As required) Steam sterilisation with an effective and validated procedure: <ol style="list-style-type: none"> Non-critical instruments[*]; steam sterilisation not required Semi-critical instruments^{**}; unpackaged on trays, sieve trays, tablets, stands, bowls, etc. or packed in sterile product packaging in accordance with EN 868 Critical instruments^{***}; packaged in sterile product packaging in accordance with EN 868 Documented release of the medical devices for use Contamination-protected storage, storage duration acc. to EN 868 <p>For general and surgical instruments:</p> <p>For rotating instruments:</p> <p>Washer-disinfector:</p> <p>Program: Temperature/hold time:</p> <p>Steam steriliser: autoclave</p> <p>Program: Temperature/hold time:</p> <p>Cabinet, drawer (protected against dust)</p>		<p>Processing is performed directly after the treatment, taking into account the risk assessment and classification of the medical devices</p> <p>dry storage in washer-disinfector before disinfection, not over 6 hrs.</p>	
<p>Transmission instruments for general, preventive, restorative, orthodontic, surgical, periodontal and endodontic treatments, e.g. turbine hand and angle pieces</p> <p>Auxiliary devices with/without leakage of fluids and/or air and particles, e.g. dental clearance equipment, powder jet equipment, intraoral camera, etc.</p>	<p>Processing is performed in the same way as for the instruments described above (steps 3 and 9 according to the manufacturer's information)</p> <p>External: wipe disinfection and cleaning</p> <p>Interior: cleaning Care (lubrication of interior parts)</p> <p>Automated reprocessing: Disinfection/cleaning e.g. in a washer-disinfector Internal cleaning and maintenance of internal parts</p> <p>Preparation, unit: Unit: Program: Temperature/hold time:</p>		<p>After every treatment and after contamination</p>	
Surfaces				
<p>Surfaces and floors in examination and treatment rooms (e.g. treatment unit, devices, cabinets)</p>	<p>Disinfect and clean surfaces with a mopping and wiping disinfection procedure (using an appropriate wet wiping method). No rinsing or wiping dry before the end of the reaction time. Re-use the surfaces when they are visibly dry.</p> <p>With single-use wipes or hygienically prepared wipes, multi-use wiping cloths</p>		<p>Every working day (e.g. at the end of the treatment day), additionally as required</p>	<p>Cleaner</p>
<p>Objects and small surfaces in the patient vicinity that have been contaminated and cannot be disinfected thermally or by immersion in solutions (e.g. treatment unit, devices, work surfaces, equipment parts, X-ray equipment, lamp handles, intraoral film shells)</p>	<p>Wiping with quick disinfectants Limit spray disinfection to areas that cannot be reached by wipe disinfection.</p> <p>Ensure complete wetting, comply with reaction time, allow to dry</p> <p>Cleaning and care (e.g. treatment chair)</p>		<p>In the treatment area after every treatment and after contamination</p> <p>As required</p>	
Skin protection plan				
<p>Skin, when affected by</p> <ul style="list-style-type: none"> water-soluble working materials, such as disinfectants, cleaning agents, X-ray development chemicals non water-soluble materials, such as organic solvents, petroleum ether accumulation of moisture and softening of skin (maceration), e.g. when wearing gloves 	<p>Cleaning hands: Under running water with gentle liquid soap from a dispenser, dry with single-use towels</p> <p>Skin protection and care: Apply carefully between fingers, paying special attention to fingernail folds and the back of the hands</p>	<p>Gentle liquid washing preparations: Single-use towel</p> <p>Protection emulsion W/O Preparation:</p>	<p>Before starting work, after finishing work, in case of visible soiling, before meals, following every visit to the toilet, after sneezing, blowing one's nose or coughing</p>	<p>All employees</p>
<p>Skin, in the event of non-physiological colonisation with e.g. bacteria, fungi, viruses</p>	<p>Hygienic hand disinfection: Hand disinfectant in accordance with the rub-in process CEN/EN 1500 to be rubbed in until dry. Hands must be air-dry before putting on the medical single-use gloves.</p> <p>Surgical hand disinfection: Wash hands and forearms up to the elbow for 1 minute, brushing fingernails and nail folds if required. Wet hands and forearms completely by rubbing in the hand disinfectant solution in portions (following the rubbing-in method described in CEN/EN 12791); keep moist during the reaction time. Hands must be air-dry before putting on sterile medical single-use gloves.</p>		<p>Before preparing for work, before and after each treatment, after interruptions in treatment, after workstation maintenance, before putting on or taking off gloves</p> <p>Before extensive dental-surgical or oral-surgical procedures, before all dental-surgical or oral-surgical procedures on patients with an increased risk of infection</p>	<p>All employees</p>
Special areas				
<p>Suction unit, suction hoses (internal)</p> <p>Spittoon (internal)</p>	<p>Cleaning by aspirating a water-air mix, especially when blood has been aspirated</p> <p>a) Disinfection: slow suction of a mix of air and disinfection/cleaning solution</p> <p>b) Cleaning: removal of deposits (calcium, bacterial contamination, salts from jet-stream equipment)</p> <p>Disinfection and cleaning of the spittoon</p> <p>Filter change or cleaning in accordance with the manufacturer's specifications</p>	<p>Cold water (at least 1/2 litre)</p> <p>Use PPE (personal protective equipment)</p>	<p>After every treatment</p> <p>a) 1 to 2 times per day, at a minimum at the end of the working day</p> <p>b) at least 2 times per week, prior to the midday break</p> <p>Additionally as required</p>	
<p>Filter</p> <p>Amalgam separator</p>	<p>Replacement/drainage of the collecting container and disposal of the separated material in accordance with the manufacturer information</p>	<p>Do not touch contaminated pieces of the separator, use PPE</p>	<p>Use PPE As required</p> <p>As required</p>	
<p>Suction hoses (external) etc.</p> <p>Spittoon (internal and external)</p>	<p>Disinfection and cleaning through wiping</p>		<p>After every treatment</p>	
<p>Impressions e.g. alginate, silicone, polyether rubber, Hydrocolloids, Polysulfide etc. as well as</p> <p>workpieces and equipment e.g. prostheses, bridges, etc., plaster models, bite registrations and bite plates made of wax, shellac or plastics, as well as articulators, occluders, fixators, face bows and more.</p>	<p>Disinfection and cleaning with Automated procedures: in accordance with the manufacturer's specifications</p> <p>Wet-chemical procedures: Disinfection and cleaning by soaking in solution, after the advised reaction time rinse under running water and air dry if necessary</p>	<p>Device/manufacturer:</p>	<p>After removal from the patient's mouth Prior to transport to the dental lab After receipt from the dental lab</p>	
<p>Objects with acid-soluble coatings e.g. cement, calculus, oxide film, flux residues</p>	<p>Cleaning (e.g. in an ultrasonic bath)</p>		<p>As required</p>	
<p>Objects with coatings of plaster or alginate</p>	<p>Cleaning (e.g. in an ultrasonic bath)</p>		<p>As required</p>	
Miscellaneous				
<p>Oral antiseptics</p>	<p>Rinse the oral cavity or soak swabs with antiseptic to reduce germs in the oral cavity of the patient while safeguarding the surgery team against infection</p>	<p>Preparation:</p>	<p>Before surgical procedures, for patients with an increased risk of infection</p>	<p>Patient</p>
<p>Practice laundry Textile practice clothing, professional clothing, towels, wiping cloths, etc.</p>	<p>Practice laundry is to be collected in sufficiently resistant and thick containers/sacks and reprocessed with thermal or chemical-thermal wash and disinfection procedures</p>	<p>Device/manufacturer: Wash program/temperature: Preparation/dosage:</p>	<p>After use After contamination</p>	<p>Operator</p>
<p>Waste Waste from treatment and examination rooms</p>	<p>Pointed, sharply, fragile objects, body parts (e.g. teeth, tissue) with blood, secretions from contaminated refuse (e.g. saliva ejector, swabs, etc.), as well as waste contaminated with agents of notifiable diseases must be immediately collected at their point of origin in non-penetrable, tear- and moisture-proof and sealed containers. The contents must not be, transferred to other containers nor sorted, and must be disposed of in securely sealed containers for burning. Lead foil, old films, fixing baths, developer, residues containing mercury, old medications and batteries are disposed of by specialist disposal firms. Household-like waste is collected in waste bags and disposed of in accordance with the local guidelines.</p>			

^{*} Non-critical medical devices: medical devices that merely come into contact with intact skin
^{**} Semi-critical medical devices: medical devices that come into contact with mucosa or pathologically affected skin
^{***} Critical medical devices: medical devices that penetrate skin or mucosa and come into contact with blood, internal tissue or organs (including wounds) in the process, and medical devices for application of blood, blood products and sterile drugs. Patients suffering from or thought to suffer from clinical conditions of transmissible spongiform encephalopathy (CJK, vCJK etc.) are to be treated in institutions that are sufficiently equipped to prevent infection. The processing of re-usable instruments and materials must only be performed in institutions with an externally certified quality management system in accordance with DIN EN ISO 13485/13488.

¹⁾ For further information, refer to the product information leaflets available on the Internet.
²⁾ In accordance with EN 13727, EN 13624, EN 14561, EN 14562
³⁾ Use disinfectants safely. Always read the label and product information before use.